

BIOSYSTEMS A25

ALKALINE PHOSPHATASE mod. IFCC

Liquid
2 Reagents

REF

Cont.

| | | | |
|---------------|-------------------|------------|-----------|
| D95564 | 5 x 100 ml | 4 x 100 ml | Reagent 1 |
| | | 1 x 100 ml | Reagent 2 |
| D95565 | 5 x 50 ml | 4 x 50 ml | Reagent 1 |
| | | 1 x 50 ml | Reagent 2 |

| | | | |
|--------|-----------|------------------|-------------|
| D98485 | 5 x 3 ml | Calibrator | Diacal Auto |
| D98481 | 12 x 5 ml | Control normal | Diacon N |
| D98482 | 12 x 5 ml | Control abnormal | Diacon P |

1. Reagent preparation:

The reagents are ready to use.

2. Instrument settings:

Temperature: 37 °C

| | | |
|-------------|----------------------------|---------------------|
| GENERAL | Test name | ALP IFCC |
| | Analysis mode | Kinetic |
| | Sample type | Serum |
| | Units | U/L |
| | Reaction type | Ascending |
| | Decimals | 0 |
| | Replicates | 1 |
| | Name of assoc. Constituent | - |
| PROCEDURE | Type of reading | Monochromatic |
| Volumes | Sample | 4 |
| | Reagent 1 | 200 |
| | Reagent 2 | 50 |
| | Washing | 1.2 |
| | Predilution factor | - |
| | Postdilution factor | 2 |
| Filters | Main | 405 |
| | Reference | |
| Times | Reading 1 | 120 |
| | Reading 2 | 420 |
| | Reagent 2 | 60 |
| CALIBRATION | Type of calibration | Multiple Calibrator |
| | Calibrator replicates | 3 |
| | Blank replicates | 1 |
| | Calibration curve | - |
| OPTIONS | Blank absorbance limit | 0,600 |
| | Kinetic blank limit | - |
| | Linearity limit | 4000 |

NOTE:

These suggested instructions and instrument parameters are to be used in conjunction with the reagent package insert and the instrument operation manual. Refer to these documents for complete instructions before performing the tests.