

## APPLICATION PROCEDURE

### KONELAB

### GOT (AST, ASAT) MODIFIED IFCC

Liquid  
2 Reagents

REF	Cont.
-----	-------

<b>D94610</b>	<b>5 x 100 ml</b>	4 x 100 ml 1 x 100 ml	Reagent 1 Reagent 2	<b>2500 tests</b> 200 µl/test
<b>D98616</b>	<b>5 x 50 ml</b>	4 x 50 ml 1 x 50 ml	Reagent 1 Reagent 2	<b>1250 tests</b> 200 µl/test
<b>DK0728*</b>	<b>5 x 50 ml</b>	4 x 50 ml 1 x 50 ml	Reagent 1 Reagent 2	<b>1250 tests</b> 200 µl/test

\* Reagent filled into Kone system bottles

additionally offered:

D98485	5 x 3 ml	Calibrator	Diacal Auto
D98481	12 x 5 ml	Control normal	Diacon N
D98482	12 x 5 ml	Control abnormal	Diacon P

#### 1. Reagent preparation

The reagents are ready to use.

#### 2. Instrument settings:

Temperature: 37 °C  
Substrate Start

<b>Test Definition:</b>			
Test type	Photometric		
Full name	AST 2		
On line name	AST 2		
Test type	Photometric		
Result unit	U/l		
Number of decimals	0		
Acceptance	AUTOMATIC		
Dilution 1 +	0		
Sample type	Serum/plasma		
Test in use	YES		
	Low	High	Units
Test Limit	0	1800	U/l
Initial Absorbance	1	2.2	A
Dilution limit	0	300	U/l
Secondary dil. 1 +	0	5	
Correction factor	1.00		
Correction bias	0.00		
<b>Calibration parameters</b>			
Calibration type	NONE		
Factor	-2000	Bias	0
Bias corr.in use	NO		
Manual QC in use	YES	Routine QC in use	NO
Acceptance	Manual	Interval	
Control	Mean	SD	
*	*	*	

<b>Test flow</b>			
Blank	NO	Antigen excess	NO
Reagent 1	AST1		
Reagent volume (µl)	160		
Disp with	WATER	Volume(µl)	0
Reagent 2	AST2		
Reagent Volume (µl)	40		
Disp with	WATER	Volume(µl)	0
Dilution with	WATER		
Sample Volume(µl)	20		
Disp with	WATER	Volume(µl)	0
Incubation Time (sec)	120		
	λ 1 (nm)	λ 2 (nm)	
	340		
Curve type	LINEAR		
Nonlinearity			
	CONC.(U/L)	20	
	Time (sec)	120	
	Point & Inter		
<b>Points and Intervals</b>	<b>3/54</b>		

\*) Data entry by the user

**NOTE:** These suggested instructions and instrument parameters are to be used in conjunction with the reagent package insert and the instrument operation manual. Refer to these documents for complete instructions before performing the tests.