

Α ΕΛΛΗΝΙΚΗ ΔΗΜΟΚΡΑΤΙΑ
ΥΠΟΥΡΓΕΙΟ ΠΑΙΔΕΙΑΣ, ΕΡΕΥΝΑΣ ΚΑΙ ΘΡΗΣΚΕΥΜΑΤΩΝ
ΙΝΣΤΙΤΟΥΤΟ ΤΕΧΝΟΛΟΓΙΑΣ ΥΠΟΛΟΓΙΣΤΩΝ ΚΑΙ ΕΚΔΟΣΕΩΝ ΔΙΑΔΙΚΤΥΟΥ
ΙΤΥΕ - ΙΤΕ

CERULOPLASMIN

Diagnostic reagent for the quantitative in vitro determination of Ceruloplasmin in human serum by turbidimetric assay

REF

Cont.

A00531 1 x 10 ml Ceruloplasmin Antibody Reagent
5 x 25 ml PEG4 Buffer

Additionally offered:

A00704 5 x 1 ml Protein Calibrator 5 level series
A00580 1 x 1 ml Protein Calibrator High
A00703 1 x 5 ml Protein Calibrator High
A00590 1 x 1 ml Protein Control
A00800 1 x 5 ml Protein Control
A08591 1 x 1 ml Protein Control Low
A08823 1 x 5 ml Protein Control Low

1. Reagent preparation

Sample: Ready to use

Reagent : Buffer, ready to use

Start Reagent: Antibody Reagent, ready to use

Calibrator: Ready to use

Controls: Ready to use

2. Instrument setting

Test	Ceruloplasmin
Unit	mg/dL
Decimals	1
Analysis type	Two point end assay 1
Main wavelength	340
Sub wave length	
Method	Turbidimetry
Correction	1.00x + 0
Calibration type	Spline
Standard	Multipoint calibration
Normal range	20 – 50
Sample	3 µL
Automatic dilution	NO
R1	300 µL
R2	45 µL
Read Main	54
Sub	31
Prozone check	NO

Dialab - www.trionmed.ru
- www.trionmed.com
www.kuban-pm.ru

Note: These suggested instructions and instrument parameters are to be used in conjunction with the reagent package insert and the instrument operation manual. Refer to these documents for complete instructions before performing tests.