

Α ΕΛΛΗΝΙΚΗ ΔΗΜΟΚΡΑΤΙΑ
ΥΠΟΥΡΓΕΙΟ ΥΓΕΙΑΣ
ΕΘΝΙΚΟ ΚΕΝΤΡΟ ΕΡΕΥΝΑΣ
ΔΙΑΓΝΩΣΤΙΚΗΣ ΚΑΙ ΑΝΤΙΣΤΡΑΤΕΓΜΑΤΟΣ

CYSTATIN C (Cys C)

2 Reagents

Diagnostic reagent for the quantitative in vitro determination of Cystatin C (Cys C) in human serum or plasma by turbidimetric assay

REF	Cont.	
906840	2 x 25 ml 1 x 10 ml	Cystatin C Buffer Cystatin C Latex Reagent
906850	1 x 25 ml 1 x 5 ml	Cystatin C Buffer Cystatin C Latex Reagent

Additionally offered:

906870	5 x 2 ml	Cystatin C Calibrator Set
906880	6 x 2 ml	Cystatin C Control Set (2 Levels)

1. REAGENT PREPARATION

The reagents are ready to use.

REAGENT1: Cystatin C Buffer

Reagent2: Cystatin C Latex Reagent

Use saline for the 0.0 calibrator.

2. INSTRUMENT SETTING

Test	Cystatin C
Unit	mg/L
Decimals	2
Analysis type	Two point end assay 1
Main wavelength	600
Sub wave length	
Method	Turbidimetry
Correction	1.00x + 0
Calibration type	Spline
Standard	Multipoint calibration
Normal range	0.4 – 1.2
Sample	3 µL
Automatic dilution	NO
R1	335 µL
R2	60 µL
Read Main	54
Sub	31
Prozone check	NO

www.trionmed.ru
www.trionmed.com
 Dialab - www.kuban-pm.ru

Note: These suggested instructions and instrument parameters are to be used in conjunction with the reagent package insert and the instrument operation manual. Refer to these documents for complete instructions before performing tests.