

# Α ΕΛΛΗΝΙΚΗ ΔΙΑΓΝΩΣΤΙΚΗ ΕΡΕΥΝΑ ΕΠΙΧΕΙΡΗΣΗ

## FERRITIN

Diagnostic reagent for the quantitative in vitro determination of Ferritin in human serum by turbidimetric immunoassay

REF	Cont.	
<b>A06551</b>	1 x 10 ml Ferritin Latex Reagent 2 x 25 ml Ferritin Buffer	200 tests/kit

Additionally offered:

A06560	5 x 1 ml Ferritin Calibrator 5 level series
A00610	1 x 1 ml Ferritin Control High
A00821	1 x 5 ml Ferritin Control High
A00570	1 x 1 ml Ferritin Control Low
A00820	1 x 5 ml Ferritin Control Low

### 1. Reagent preparation:

Sample:	Ready for use
Buffer:	Ready for use
Antibody reagent:	Ready for use
Calibrator:	Dilute the Protein calibrator high successively 1:2 in saline 9 g/L to set up a calibration curve. Alternatively use the ready for use Protein Calibrator Set. Use Saline 9 g/L as zero point.

### 2. Instrument setting

Test	Ferritin
Unit	ng/mL
Decimals	1
Analysis type	Two point end assay 1
Main wavelength	600
Sub wave length	
Method	Turbidimetry
Correction	1.00x + 0
Calibration type	Spline
Standard	Multipoint calibration
Normal range	10 – 300
Sample	40 µL
Automatic dilution	NO
R1	220 µL
R2	50 µL
Read Main	54
Sub	31
Prozone check	NO

:

-

[www.trionmed.ru](http://www.trionmed.ru)  
 - [www.trionmed.com](http://www.trionmed.com)  
 Dialab - [www.kuban-pm.ru](http://www.kuban-pm.ru)

**Note:** These suggested instructions and instrument parameters are to be used in conjunction with the reagent package insert and the instrument operation manual. Refer to these documents for complete instructions before performing tests.