

Α ΕΛΛΗΝΙΚΗ ΔΗΜΟΚΡΑΤΙΑ  
ΥΠΟΥΡΓΕΙΟ ΥΓΕΙΑΣ  
ΕΘΝΙΚΟ ΚΕΝΤΡΟ ΕΡΕΥΝΑΣ ΚΑΙ ΔΙΑΓΝΩΣΤΙΚΗΣ ΒΙΟΪΑΤΡΙΚΗΣ

## FIBRONECTIN

Diagnostic reagent for the quantitative in vitro determination of Fibronectin in human plasma by turbidimetric assay

REF

Cont.

A00518 1 x 5 ml Fibronectin Antibody Reagent  
2 x 25 ml PEG6 Buffer

Additionally offered:

A02712 1 x 1 ml Fibronectin Calibrator

### 1. Reagent preparation

Sample: Ready to use (plasma)

Reagent : Buffer, ready to use

Start Reagent: Antibody Reagent, ready to use

Calibrator: Ready to use

Controls: Ready to use

### 2. Instrument setting

Test	Fibronectin
Unit	mg/dL
Decimals	1
Analysis type	Two point end assay 1
Main wavelength	340
Sub wave length	
Method	Turbidimetry
Correction	1.00x + 0
Calibration type	Spline
Standard	Multipoint calibration
Normal range	25 – 40

Sample	7 µL
Automatic dilution	NO
R1	200 µL
R2	40 µL
Read Main	54
Sub	31
Prozone check	NO

\* Standard value

\*\* Data entered by operator

125 Tests/Kit

:

-

Dialab -

[www.trionmed.ru](http://www.trionmed.ru)  
- [www.trionmed.com](http://www.trionmed.com)  
[www.kuban-pm.ru](http://www.kuban-pm.ru)

**Note:** These suggested instructions and instrument parameters are to be used in conjunction with the reagent package insert and the instrument operation manual. Refer to these documents for complete instructions before performing tests.