

Α ΕΛΛΗΝΙΚΗ ΔΗΜΟΚΡΑΤΙΑ
ΥΠΟΥΡΓΕΙΟ ΥΓΕΙΑΣ ΚΑΙ ΚΟΙΝΩΝΙΚΗΣ ΑΣΦΑΛΙΣΗΣ

LIPOPROTEIN (a)

Diagnostic reagent for the quantitative in vitro determination of Lp(a) in human serum by turbidimetric assay

REF

A00527	1 x 10 ml Lp(a) Latex Reagent 5 x 25 ml Lp(a) Buffer
A00733	1 x 1 ml Lp(a) Calibrator High
A00734	1 x 5 ml Lp(a) Calibrator High
A03736	4 x 1 ml Lp(a) Calibrator 4 level series
A00816	1 x 1 ml Lp(a) Control High
A00817	1 x 5 ml Lp(a) Control High
A00814	1 x 1 ml Lp(a) Control Low
A00815	1 x 5 ml Lp(a) Control Low

1. Reagent preparation

Sample:	Ready for use
R1:	Ready for use
R2:	Ready for use
Calibration:	Dilute the Lp(a) Calibrator high successively 1:2 in NaCl 9 g/L to set up a calibration curve. Use NaCl 9 g/L as zero point. Alternatively use the ready for use Lp(a) Calibrator set.

2. Instrument setting

Test	Lipoprotein (a)
------	-----------------

Unit	mg/dL
Decimals	1
Analysis type	Two point end assay 1
Main wavelength	600
Sub wave length	
Method	Turbidimetry
Correction	1.00x + 0
Calibration type	Spline
Standard	Multipoint calibration
Normal range	0 – 30

Sample	12 µL
Automatic dilution	1:10
R1	225 µL
R2	20 µL
Read Main	54
Sub	31
Prozone check	NO

500 Tests/Kit

:

-

www.trionmed.ru
- www.trionmed.com
Dialab - www.kuban-pm.ru

Note: These suggested instructions and instrument parameters are to be used in conjunction with the reagent package insert and the instrument operation manual. Refer to these documents for complete instructions before performing tests.