

# Α ΕΛΑΤΙΟΕΥΡΩΠΟΙΟΑ ΕΕΑΥÜÒÙΝΟΘÒΑΓ Δ

## MICROALBUMIN

Diagnostic reagent for the quantitative in vitro determination of Microalbumin in human urine by turbidimetric assay

REF
-----

Cont.
-------

A00528      1 x 10 ml Microalbumin Antibody Reagent      625 tests/kit  
 5 x 25 ml Microalbumin Buffer

Additionally offered:

A00955      5 x 1 ml Microalbumin Calibrator 5 level series  
 A00707      1 x 1 ml Microalbumin Calibrator  
 A00708      1 x 5 ml Microalbumin Calibrator  
 A00808      1 x 1 ml Microalbumin Control  
 A00809      1 x 5 ml Microalbumin Control

### 1. Reagent preparation

Sample:            Ready for use  
 R1:                Ready for use  
 R2:                Ready for use  
 Calibration:      Dilute the Microalbumin calibrator successively 1:2 in NaCl 9 g/L to set up a calibration curve. Use NaCl 9 g/L as zero point. Alternatively use the Microalbumin calibrator set.

### 2. Instrument setting

Test	Microalbumin
Unit	mg/L
Decimals	1
Analysis type	Two point end assay 1
Main wavelength	340
Sub wave length	
Method	Turbidimetry
Correction	1.00x + 0
Calibration type	Spline
Standard	Multipoint calibration
Normal range	0 – 25
Sample	6 µL
Automatic dilution	NO
R1	200 µL
R2	35 µL
Read Main	[52] - [54]
Sub	[28] - [30]
Prozone check	NO

:

-

Dialab - [www.trionmed.ru](http://www.trionmed.ru)  
[www.trionmed.com](http://www.trionmed.com)  
[www.kuban-pm.ru](http://www.kuban-pm.ru)

**Note:** These suggested instructions and instrument parameters are to be used in conjunction with the reagent package insert and the instrument operation manual. Refer to these documents for complete instructions before performing tests.