

Biosystems A25

PHOSPHORUS INORGANIC MOLYBDATE

Liquid
Single Reagent

REF

Cont.

D00362 **5 x 50 ml** Single Reagent
D95362 1 x 3 ml Phosphorus Standard

D98485 5 x 3 ml Calibrator DiaCAL Auto
D98481 12 x 5 ml Control normal Diacon N
D98482 12 x 5 ml Control abnormal Diacon P

1. Reagent preparation:

The reagent is ready to use.

2. Instrument settings:

Temperature: 37 °C

GENERAL	Test name	Phosphorus
	Analysis mode	Endpoint
	Sample type	Serum
	Units	mg/dl
	Reaction type	Ascending
	Decimals	1
	Replicates	1
	Name of assoc. Constituent	-
PROCEDURE	Type of reading	Monochromatic
Volumes	Sample	3
	Reagent 1	300
	Reagent 2	
	Washing	1.2
	Predilution factor	-
	Postdilution factor	2
Filters	Main	340
	Reference	
Times	Reading 1	300
	Reading 2	
	Reagent 2	
CALIBRATION	Type of calibration	Multiple Calibrator
	Calibrator replicates	3
	Blank replicates	1
	Calibration curve	-
OPTIONS	Blank absorbance limit	
	Kinetic blank limit	-
	Linearity limit	15

NOTE:

These suggested instructions and instrument parameters are to be used in conjunction with the reagent package insert and the instrument operation manual. Refer to these documents for complete instructions before performing the tests.

:

-

Dialab - www.trionmed.ru
- www.trionmed.com
www.kuban-pm.ru